Analysis of O₂ Column Abundances using A-Band High Resolution Solar Absorption Spectroscopy at Table Mountain, CA

Z. Yang¹, P. Wennberg¹, G. Toon², L. Brown², C. Miller³, T. Pongetti², R. Cageao² and <u>S. Sander^{2,1}</u>


¹Division of Geological and Planetary Science California Institute of Technology Pasadena, CA 91125

²Jet Propulsion Laboratory California Institute of Technology Pasadena, CA 91109

³ Haverford College, Philadelphia, PA


BACKGROUND AND OBJECTIVES

- ESSP/Orbiting Carbon Observatory (OCO) measurement requirement is CO₂ column abundance with 0.3% precision.
- Ground-based FT spectroscopy plays an important role in calibration and validation of OCO.
- FTUVS measures O₂ and CO₂ in the same near-IR bands as OCO. Diurnal measurements will establish the feasibility of achieving 0.3% precision.
- O₂ is the first objective because its column abundance is known.
- Provide a check on A-band spectroscopic parameters that will be useful for other instruments including SAGE III, SCIAMACHY




OCO IFOV and spectral bands

Table mountain dome and instruments







Atmospheric Measurements


- O₂ A-band data acquired 3/24/03 4/11/03 using FTUVS/TMF (12days)
- Air mass range 7.7 (morning) 7.6 (afternoon)
- Spectra every 30 min (A/M < 3), continuously (A/M > 3)
- Spectral range: 12,100 15100 cm⁻¹, resolution 0.056 cm⁻¹
- Noontime SNR: 1300:1


Forward Model Configuration

- O₂ line positions and strengths from Brown and Plymate (2000)
- Update B&P's broadening coefficients with high N' results of Yang et al.
- Super-Lorentz lineshape of Hirono and Nakazawa (1982) to model continuum
- Temperature profile:
 - 2.3 < z < 5 km: best fit to spectrum
 - z > 5 km: profile from assimilation model (NCEP/NCAR reanalysis)
- Atmospheric overburden calculated from measured surface pressure and assumption of hydrostatic equilibrium
- Water vapor is not explicitly retrieved; climatological water vapor is used as a first-order correction for the column-average O₂ mixing ratio.
- O₂ isotopomers are retrieved separately.


Lorentz:
$$k(v) = \frac{S}{\pi} \frac{\gamma}{(v - v_0)^2 + \gamma^2}$$

Super-Lorentz:
$$k(\nu) = \frac{S}{\pi^{1/2}} \frac{S\Gamma(\eta/2)}{\Gamma((\eta-1)/2)} \frac{\gamma^{\eta-1}}{\{(\nu-\nu_0)^2 + \gamma^2\}^{\eta/2}}, \quad (1 < \eta \le 2)$$







O₂ A-Band Spectrum – TMF/FTUVS



TMF/FTUVS O_2 A-band: 3/24/03 - 4/10/03 super-Lorentz line shape, eta = 1.977



TMF/FTUVS O_2 A-band: 3/24/03 - 4/10/03 super-Lorentz line shape, eta = 1.958




Conclusions

• The O_2 column abundance can be measured with a precision of $\pm 0.3\%$ using ground-based high resolution solar absorption spectroscopy in the A-band.

First retrieval of an atmospheric gas with 0.3% measurement precision.

- Serious spectroscopic issues remain:
 - * Characterization of line mixing
 - * Line shape in the far wings.
 - * Line strengths, air- and self-broadened half-widths esp. at high J.
 - * Spectroscopic parameters for isotopomers.
- Other retrieval issues:
 - * Temperature profile.
 - * Effects of clouds; multiple scattering?


Air and Self-Broadened Half-Widths vs. N' for the O₂ A-Band



FTUVS INSTRUMENT INSTALLATION AT TABLE MOUNTAIN FACILITY



View of the building which houses the FTUVS instrument at JPL's Table Mountain Facility at 7500 ft. elevation near Wrightwood, California.



Layout of the FTUVS facility showing the tracker (heliostat), telescope and FTUVS interferometer.



The FTUVS heliostat showing the 20 in. primary mirror and the 12.6 in. secondary mirror. The secondary can be easily removed from its mounting for zenith sky observations.



View of the FTUVS interferometer (inside the large white box) and the telescope. The telescope achieves a 3:1 magnification of the solar image and has a CCD camera for fine tracking which views the solar image reflected off the field stop.